# DEI® Rainbow

---

# Hydro Big One Phase

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**DEI® Rainbow**

# Hydro Big One Phase

**MONOPHASE COMPOUND  
FOR UPMOST HYDROCOMPATIBLE  
PRECISION IMPRESSIONS**

Product for dental use  
Keep away from children  
Date of the specifications: 10/2009

#### Available packages

Two 380 ml cartridges + 10 cannulas (**COD. RAINH 10**)

#### Characteristics

DEI® Rainbow Hydro Big One Phase is a high viscosity material based on vinylpolisiloxane technique, indicated for bridges, crowns, inlays, relinings and in technique of copper caps impression.

DEI® Rainbow Hydro Big One Phase is packed in cartridges for mixing devices.

It conforms to ISO 4823 type 1 regulation.

Shore A about 65.

#### Advantages

DEI® Rainbow Hydro Big One Phase shows high dimensional stability, hydrocompatibility and elastic memory, is odourless, tasteless, has an high tear resistance and is thixotropic (stable on impression tray, however flowing when undergone to pressure).

#### Instructions

Extrude the product on the impression tray. Insert the impression tray into the oral cavity within 1 minute from the extrusion, and let it in place at least 2 minutes.

The indicated times refer to room temperature of 23°C.

Higher temperature shorten the times, lower ones lengthen them.

#### Warnings

The impression can be cast 30 minutes after the taking away from oral cavity.

It can be electroplated in copper or silver.

Disinfect with suitable solutions. DEI® Italia Sicursept is recommended.

Use specific cannulas for every type of cartridge.

#### Preparation of the cartridge

Take away the safety cap, insert the cannula.

Block turning the supplied nut half a turn, insert the cartridge into the mixing device.

Don't use again the safety cap.

After use leave the cannula on the cartridge.

Only for new cartridges: for better safety, discard the first 2-3 cm of dispensed compound.

Via Torino 765 - 21020 Mercallo (VA)  
Phone +39 0331 969270 - Fax +39 0331 969271  
www.deiitalia.it

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